



**Technical Data Sheet**

**Model:** DR.LNS1717 & DR.LNS1717WD

**Description:** Panel Protector

**Product Attributes:**

- Lightweight
- Available with a grid or carbon fiber panel for
- Multiple handle orientations: Long, Short, and Dual
- Multiple layers of impact protection
- For cassette-sized wireless or tethered panels

**Technical Specifications**

Model	DR.LNS1717 & DR.LNS1717WD
Weight	2.08 lbs [0.94 kg] Plus Grid 3.58 lbs [1.62 kg] Without Grid
Size	21.84" X 19.38" X .97" (554.7mm X 492.2 mm X 24.5mm)
Box Dimensions	24" X 24" X 7" (609.6mm X 609.6mm X 177.8mm)
Shipping Weight	9 lbs. (4.1 kg)

**Statement of Conformance**

RCI Operating CO LLC d/b/a RC Imaging hereby states the data, specifications and descriptions contained herein this document are accurate and representative of the model described.

**Registration and Compliance**

The models listed above are made in the USA and registered to:

RCI Operating CO LLC  
50 Old Hojack Lane  
Hilton, NY 14468

**Device Family Classification**

COUNTRY	CLASS
USA	FDA, Class 1, 510 (k) exempt
Canada	Class 1
Europe	Class 1
UK	Class 1