



### Technical Data Sheet

**Model:** WB.C1012

**Description:** Weight Bearing Cover

#### Product Attributes:

- Light Weight
- High-density polyethylene frame
- Polycarbonate Top
- Easy lift handle
- Rubber Feet

#### Technical Specifications

Model	WB.C1012
Weight	4.3 lbs. (1.9 kg)
Size	16.43" X 12.90" X 1.52" (417.4mm X 327.7mm X 38.5mm)
Box Dimensions	20" X 17" X 6" (508mm X 431.8mm X 152.4mm)
Shipping Weight	6 lbs. (2.8 kg)

#### Statement of Conformance

RCI Operating CO LLC d/b/a RC Imaging hereby states the data, specifications and descriptions contained herein this document are accurate and representative of the model described.

#### Registration and Compliance

The model listed above is made in the USA and is registered to:

RCI Operating CO LLC  
50 Old Hojack Lane  
Hilton, NY 14468

#### Device Family Classification

COUNTRY	CLASS
USA	FDA, Class 1, 510 (k) exempt
Canada	Class 1
Europe	Class 1
UK	Class 1