

**Technical Data Sheet**

**Model:** WB.C1417

**Description:** Weight Bearing Cover

**Product Attributes:**

- Light Weight
- High-density polyethylene frame
- Smoked Polycarbonate top
- Easy lift handle
- Rubber Feet

**Technical Specifications**

<b>Model</b>	<b>WB.C1417</b>
Weight	7.14 lbs. (3.24 kg)
Size	21.43" X 16.90" X 1.52" (544.4mm X 429.3mm X 38.5mm)
Box Dimensions	23" X 20" X 6" (584.2mm X 508mm X 152.4mm)
Shipping Weight	11 lbs. (5.0 kg)

**Weight Tolerance**

Weight Tolerance: 500 lbs. (226.80 kg)

**Statement of Conformance**

RCI Operating CO LLC d/b/a RC Imaging hereby states the data, specifications and descriptions contained herein this document are accurate and representative of the model described.

**Registration and Compliance**

The model listed above is made in the USA and is registered to:

RCI Operating CO LLC  
50 Old Hojack Lane  
Hilton, NY 14468

**Device Family Classification**

<b>COUNTRY</b>	<b>CLASS</b>
USA	FDA, Class 1, 510 (k) exempt
Canada	Class 1
Europe	Class 1
UK	Class 1