

Weight Bearing Tunnel



This product is ideal for equine and veterinary imaging to protect your DR Panel from excessive weight.

Weight Bearing Tunnel

Our Weight Bearing Tunnel is perfect for imaging large animals when paired with one of our patented Lock-N-Secure® panel protectors, this product will allow for seamless removal without disruption.



Great Features

- Distributed weight capacity ranging from 700lbs to over 2000lbs
- Grab-N-Go handle
- Target area for easy alignment

Quality Construction

- Available in Carbon Fiber and Polycarbonate
- Non-skid rubber feet

Available in Multiple Sizes

- 10" x 12" (25 x 30cm)
- 14" x 17" (35 x 43cm)
- 24 x 30cm



RC Imaging
your radiology resource

800.528.6311 585.392.4336 585.392.5575
50 Old Hojack Lane, Hilton, New York 14468
www.rcimaging.com

