

**Product Datasheet**

**Model:** DR.LNS1417SD2 & DR.LNS1417SD2WD


**Description:** Panel Protector

**Product Attributes:**

- Lightweight
- Available with a grid or carbon fiber panel for
- Multiple handle orientations: Long, Short, and Dual
- Multiple layers of impact protection
- For cassette-sized wireless or tethered panels

**Technical Specifications**

Model	DR.LNS1417SD2 & DR.LNS1417SD2WD
Weight	2.16 lbs [0.97 kg] Plus Grid 3.3 lbs [1.49 kg] Without Grid
Size	21.84" X 17.63" X .97" (554.7mm X 447.6 mm X 24.5mm)
Box Dimensions	23" X 20" X 6" (584.2mm X 508mm X 152.4mm)
Shipping Weight	9 lbs. (4.1 kg)



P/N: DR.LNS1417SD2WD      Weight: 3.3 lbs (1.49 kg)  
 P/N: DR.LNS1417SD2      Weight: 2.16 lbs (.97 kg) Plus Grid

**Statement of Conformance**

RC Imaging hereby states the data, specifications and descriptions contained herein this document are accurate and representative of the model described.

**Registration and Compliance**

The models listed above are made in the USA and registered to:

RC Imaging  
50 Old Hojack Lane  
Hilton, NY 14468

**Device Family Classification**

<b>COUNTRY</b>	<b>CLASS</b>
USA	FDA, Class 1, 510 (k) exempt
Canada	Class 1
Europe	Class 1
UK	Class 1