

Technical Data Sheet

Model: A.LPS

Description: Leg Positioning Stand

Product Attributes:

- Cam-Stay technology with 10 locking positions.
- Closed cell PSA foam for easy cleaning and sanitation.
- Non-skid rubber feet for table grid.
- Wide powder coated steel base for patient stability.
- Height range of 13.47 in. – 19.00 in.
- Light Weight

Technical Specifications

Model	A.LPS
Weight	5.6 Lbs. (2.5 Kg.)
Size	10 x 10 x 14-19 in. (254 mm. x 254 mm. x 342-483 mm)
Box Dimensions	11 x 11 x 16 ½ in. (279 x 279 x 419 mm)
Shipping Weight	8 Lbs. (3.63 Kg.)

Statement of Conformance

RCI Operating CO LLC d/b/a RC Imaging hereby states the data, specifications and descriptions contained herein this document are accurate and representative of the model described.

Registration and Compliance

The model listed above is made in the USA and is registered to:

RCI Operating CO LLC
50 Old Hojack Lane
Hilton, NY 14468

Device Family Classification

COUNTRY	CLASS
USA	FDA, Class 1, 510 (k) exempt
Canada	Class 1
Europe	Class 1
UK	Class 1