



# RC Imaging

your radiology resource

## Leg Positioning Stand



*RC Imaging Leg Positioning Stand converts to 10 telescopic locking positions for effortless vertical adjustments. Lined with pressure sensitive closed cell foam for patient comfort, non-skid rubber feet, and wide steel base designed for added stability. Ideal for obtaining Dunn view images and upright PCL stress views.*

- *Cam-Stay technology with 10 locking positions.*
- *Closed cell PSA foam for easy cleaning and sanitation.*
- *Non-skid rubber feet for table grip.*
- *Wide powder coated steel base for patient stability.*
- *Weighs 5.6 lbs.*
- *Base Plate: 10 x 10 x .25 in.*
- *Min Height: 13.47 in.*
- *Max Height: 19.00 in.*
- *Leg Holder: 7.5 x 6 x .25 in.*