

**Technical Data Sheet**

**Model:** A.MPP

**Description:** Mobile Positioner Pro DR panel holder accommodates DR panels from 10x12 in. to 17x17 in., 1.25 in. channel depth.

**Product Attributes:**

- Fully adjustable height from floor level up to 73 in. with 26 locking positions
- Rotates 360 degrees, tilts +/- 90 degrees allowing cassettes/DR panels to face straight up, down, or anywhere in between
- Easy-glide vertical and horizontal adjustable clamps extend to 23 in. in all directions
- Compatible with DR panels and cassettes from 8x10 in. to 17x17 in. including DR Panel Protectors and Grid Encasements
- 1.25 in. Channel depth
- 4 locking casters for secure positioning
- Optional 40-pound counterweighted base for increased stability and extension during cross-table procedures
- Easy assembly.

**Technical Specifications**

<b>Model</b>	<b>A.MPP</b>
<b>Weight</b>	45 Lbs. (22 kg)
<b>Size</b>	21.4x25.3x66 in. (543.2x643.22x1677mm)
<b>Box Dimensions</b>	Box 1- 70x2x2 in. (1778x50.8x50.8 mm) Box 2- 30x25x10 in. (762x635x254 mm)
<b>Shipping Weight</b>	Box 1- 9 Lbs. (4.08 kg) Box 2- 38 Lbs. (17.24 kg)

**Statement of Conformance**

RCI Operating CO LLC d/b/a RC Imaging hereby states the data, specifications and descriptions contained herein this document are accurate and representative of the model described.

**Registration and Compliance**

The model listed above is made in the USA and is registered to:

RCI Operating CO LLC  
50 Old Hojack Lane  
Hilton, NY 14468

**Device Family Classification**

<b>COUNTRY</b>	<b>CLASS</b>
USA	FDA, Class 1, 510 (k) exempt
Canada	Class 1
Europe	Class 1
UK	Class 1