

Technical Data Sheet

Model: A.MPPCW

Description: Mobile Positioner Pro DR panel holder accommodates DR panels from 10x12 in. to 17x17 in. with counter weight for extended reach, 1.25 in. channel depth.

Product Attributes:

- Fully adjustable height from floor level up to 73 in. with 26 locking positions
- Rotates 360 degrees, tilts +/- 90 degrees allowing cassettes/DR panels to face straight up, down, or anywhere in between
- Easy-glide vertical and horizontal adjustable clamps extend to 23 in. in all directions
- Compatible with DR panels and cassettes from 8x10 in. to 17x17 in. including DR Panel Protectors and Grid Encasements
- 1.25 in. Channel depth
- 4 locking casters for secure positioning
- Optional 40-pound counterweighted base for increased stability and extension during cross-table procedures
- Easy assembly.

Technical Specifications

Model	A.MPPCW
Weight	83 Lbs. (38 kg)
Size	22.2x25.3x66 in. (564x63.2x1677 mm)
Box Dimensions	Box 1- 70x2x2 in. (1778x50.8x50.8 mm) Box 2- 30x25x10 in. (762x635x254 mm) Box 3- 20x6x6 in. (508x152.4x152.4 mm)
Shipping Weight	Box 1- 9 Lbs. (4.08 kg) Box 2- 38 Lbs. (17.24) Box 3- 40 Lbs. (18.14 kg)

Statement of Conformance

RCI Operating CO LLC d/b/a RC Imaging hereby states the data, specifications and descriptions contained herein this document are accurate and representative of the model described.

Registration and Compliance

The model listed above is made in the USA and is registered to:

RCI Operating CO LLC
50 Old Hojack Lane
Hilton, NY 14468

Device Family Classification

COUNTRY	CLASS
USA	FDA, Class 1, 510 (k) exempt
Canada	Class 1
Europe	Class 1
UK	Class 1